

RPES8081

Product Information

Product SKU:	RPES8081	Expression Host:	E.coli	Size:	20µg
Tag:	N-Trx	Reactivity:	Human	Accession:	Q8NAU1

Additional Information

Calculated MW:	23.5 kDa	Observed MW:	25 kDa
Sequence:	Met144-lie212		

Protein Information

Background: -

Synonyms: FRCP, FNDC, Fibronectin type III domain-containing protein, Fibronectin type III repeat-containing protein, FRCP2, FNDC5, Irisin, Fibronectin type III domain-containing protein 5, Fibronectin type III repeat-containing protein 2

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.