

**RPES8143**

---

## Product Information

<b>Product SKU:</b>	RPES8143	<b>Expression Host:</b>	E.coli	<b>Size:</b>	20µg
<b>Tag:</b>	N-Sumo	<b>Reactivity:</b>	Human	<b>Accession:</b>	P11308-4

---

## Additional Information

<b>Calculated MW:</b>	50.2 kDa	<b>Observed MW:</b>	46 kDa
<b>Sequence:</b>	Aeg140-Tyr479		

---

## Protein Information

**Background:** -

**Synonyms:** -

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Purity:** > 90% as determined by reducing SDS-PAGE.

**Bio-Activity:** Not validated for activity

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.