

RPES8146

Product Information

Product SKU:	RPES8146	Expression Host:	E.coli	Size:	20µg
Tag:	N-Trx	Reactivity:	Mouse	Accession:	P84444

Additional Information

Calculated MW:	32 kDa	Observed MW:	35 kDa
Sequence:	Ser24-Gly133		

Protein Information

Background: -

Synonyms:	6Ckine, Beta-chemokine exodus-2, C-C motif chemokine 21a, CC21A_MOUSE, Ccl21a, Small-inducible cytokine A21a, TCA4, Thymus-derived chemotactic agent 4
Endotoxin:	< 10 EU/mg of the protein as determined by the LAL method
Formulation:	Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.
Purity:	> 90% as determined by reducing SDS-PAGE.
Bio-Activity:	Not validated for activity
Storage:	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.