

RPES8179

Product Information

Product SKU:	RPES8179	Expression Host:	E.coli	Size:	20µg
Tag:	N-Sumo	Reactivity:	Mouse	Accession:	O35228

Additional Information

Calculated MW:	36 kDa	Observed MW:	38 kDa
Sequence:	Tyr19-Pro228		

Protein Information

Background: -

Synonyms: EBI, Lymphopoietin 1(LP-1), pre-B-cell factor, Il27b, Ebi3

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.