

RPES8189

Product Information

Product SKU:	RPES8189	Expression Host:	E.coli	Size:	20µg
Tag:	N-Trx	Reactivity:	Mouse	Accession:	Q9WV56

Additional Information

Calculated MW:	32 kDa	Observed MW:	34 kDa
Sequence:	Ser319-Arg428		

Protein Information

Background: -

Synonyms: Bmp, GDF, Bmp9, Gdf2, BMP-9, Bone Morphogenetic Protein 9

Endotoxin: < 10 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.