

RPES8193

---

**Product Information**

<b>Product SKU:</b>	RPES8193	<b>Expression Host:</b>	E.coli	<b>Size:</b>	20µg
<b>Tag:</b>	N-Trx	<b>Reactivity:</b>	Mouse	<b>Accession:</b>	Q9Z126

---

**Additional Information**

<b>Calculated MW:</b>	28.2 kDa	<b>Observed MW:</b>	34 kDa
<b>Sequence:</b>	Val30-Ser105		

---

**Protein Information**

**Background:** -

**Synonyms:** Scyb, Cxcl, Cxcl4, Scyb4, Pf4

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Purity:** > 90% as determined by reducing SDS-PAGE.

**Bio-Activity:** Not validated for activity

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.