

**RPES8201**

---

## Product Information

<b>Product SKU:</b> RPES8201	<b>Expression Host:</b> E.coli	<b>Size:</b> 20µg
<b>Tag:</b> N-Trx	<b>Reactivity:</b> Mouse	<b>Accession:</b> P35461

---

## Additional Information

<b>Calculated MW:</b> 28.5 kDa	<b>Observed MW:</b> 31 kDa
<b>Sequence:</b> Leu27-Asn105	

---

## Protein Information

**Background:** -

**Synonyms:** Gr1, Gr-1, Ly-6G, Ly-6G.1, Lymphocyte antigen 6 complex locus G, Lymphocyte antigen 6G, Ly6g

**Endotoxin:** < 10 EU/mg of the protein as determined by the LAL method

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Purity:** > 90% as determined by reducing SDS-PAGE.

**Bio-Activity:** Not validated for activity

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.