

RPES8274

Product Information

Product SKU: RPES8274	Expression Host: Mammalian	Size: 20µg
Tag: C-Fc	Reactivity: Human	Accession: Q14956-2

Additional Information

Calculated MW: 75.9 kDa	Observed MW: 80-100 kDa
Sequence: Lys23-Asn486	

Protein Information

Background: -

Synonyms: -

Endotoxin: < 1.0 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.