

Recombinant Human Follistatin-related Gene Protein/FLRG Protein (His Tag)

RPES8389



Product Information

Product SKU: RPES8389	Expression Host: Mammalian	Size: 20µg
Tag: C-His	Reactivity: Human	Accession: O95633

Additional Information

Calculated MW: 25.9 kDa	Observed MW: 38 kDa
Sequence: Met27-Val263	

Protein Information

Background: -

Synonyms: UNQ, PRO, FSTL, FLRG, FSTL3, UNQ674, PRO1308

Endotoxin: < 1.0 EU/mg of the protein as determined by the LAL method

Formulation: Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

Purity: > 90% as determined by reducing SDS-PAGE.

Bio-Activity: Not validated for activity

Storage: Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.