

**RPES8408**

---

## Product Information

<b>Product SKU:</b>	RPES8408	<b>Expression Host:</b>	Mammalian	<b>Size:</b>	20µg
<b>Tag:</b>	C-Fc	<b>Reactivity:</b>	Mouse	<b>Accession:</b>	Q9JIP3

---

## Additional Information

<b>Calculated MW:</b>	54.4 kDa	<b>Observed MW:</b>	60-80 kDa
<b>Sequence:</b>	Arg18-Gly286		

---

## Protein Information

**Background:** -

**Synonyms:** -

**Endotoxin:** < 1.0 EU/mg of the protein as determined by the LAL method

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Purity:** > 90% as determined by reducing SDS-PAGE.

**Bio-Activity:** Not validated for activity

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.