

**RPES8439**

---

## Product Information

<b>Product SKU:</b>	RPES8439	<b>Expression Host:</b>	Mammalian	<b>Size:</b>	20µg
<b>Tag:</b>	C-Fc	<b>Reactivity:</b>	Mouse	<b>Accession:</b>	Q9QYY7

---

## Additional Information

<b>Calculated MW:</b>	40.8 kDa	<b>Observed MW:</b>	45 kDa
<b>Sequence:</b>	Ala22-Arg184		

---

## Protein Information

**Background:** ESM1 is a secreted protein which is produced by adipocytes. It has been noticed that ESM1 may play some role in obesity-associated vascular disease since circulating ESM-1 levels are reduced in the overweight and obese. ESM1 is mainly expressed in the endothelial cells in human lung and kidney tissues. The expression of ESM1 gene is regulated by cytokines, suggesting that it may play a role in endothelium-dependent pathological disorders. Recently, ESM1 has been described as a specific biomarker of tip cells during neoangiogenesis. Its expression has been shown to be increase in presence of pro-angiogenic growth factors such as VEGF or FGF-2. In hypervascularized cancers, overexpression of endocan has been detected by immunohistochemistry using monoclonal antibodies against ESM1.

**Synonyms:** endocan, Endothelial cell-specific molecule 1 ,

**Endotoxin:** < 1.0 EU/mg of the protein as determined by the LAL method

**Formulation:** Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.

**Purity:** > 90% as determined by reducing SDS-PAGE.

**Bio-Activity:** Not validated for activity

**Storage:** Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.