

RPES8471

Product Information

Product SKU:	RPES8471	Expression Host:	Mammalian	Size:	20µg
Tag:	C-His	Reactivity:	Mouse	Accession:	Q9JI71

Additional Information

Calculated MW:	54.7 kDa	Observed MW:	80 kDa
Sequence:	Ser28-Pro525		

Protein Information

Background: Delta-like protein 4 (DLL4 , Delta4) , a type I membrane-bound Notch ligand , is one of five known Notch ligands in mammals and interacts predominantly with Notch 1 , which has a key role in vascular development. Recent studies yield substantial insights into the role of DLL4 in angiogenesis. DLL4 is induced by vascular endothelial growth factor (VEGF) and acts downstream of VEGF as a 'brake' on VEGF-induced vessel growth , forming an autoregulatory negative feedback loop inactivating VEGF. DLL4 is downstream of VEGF signaling and its activation triggers a negative feedback that restrains the effects of VEGF. Attenuation of DLL4/Notch signaling results in chaotic vascular network with excessive branching and sprouting. DLL4 is widely distributed in tissues other than vessels including many malignancies. Furthermore , the molecule is internalized on binding its receptor and often transported to the nucleus. In pathological conditions , such as cancer , DLL4 is up-regulated strongly in the tumour vasculature. Blockade of DLL4-mediated Notch signaling strikingly increases nonproductive angiogenesis , but significantly inhibits tumor growth in preclinical mouse models. In preclinical studies , blocking of DLL4/Notch signaling is associated with a paradoxical increase in tumor vessel density , yet causes marked growth inhibition due to functionally defective vasculature. Thus , DLL4 blockade holds promise as an additional strategy for angiogenesis-based cancer therapy.

Synonyms: DLL, DII4, Delta4

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.

Endotoxin:	< 1.0 EU/mg of the protein as determined by the LAL method
Formulation:	Lyophilized from a 0.2 µm filtered solution in PBS with 5% Trehalose and 5% Mannitol.
Purity:	> 90% as determined by reducing SDS-PAGE.
Bio-Activity:	Not validated for activity
Storage:	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.