

HDBS0167

Product Information

Product SKU:	HDBS0167	Host Species:	Chimeric	Size:	100 µg
Target:	FZD10	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	Q9ULW2
Isotype:	IgG1(E239D, M241L)		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD350, FZ-10, Fz10, FzE7, hFz10

Common Name: OTSA101, ONC-005

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.