

HDBS0183

Product Information

Product SKU: HDBS0183 **Host Species:** Chimeric Rabbit **Size:** 100 µg
Target: PD-1 **Reactivity:** Canine

Additional Information

Conjugate: Unconjugated **Uniprot:** NP_001301026.1
Isotype: IgG

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1, PDCD1, CD279, SLEB2, hPD-1, hPD-I, hSLE1

Common Name: 4F12

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.