

HDBS0193

Product Information

Product SKU: HDBS0193 **Host Species:** Homo sapiens **Size:** 100 µg
Target: ANGPTL3 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** Q9Y5C1
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: ANL3; ANG-5; FHBL2; ANGPT5

Common Name: REGN 1500, evinacumab-dgnb

Recommended Dilutions: ELISA 1:5000-10000

Applications: ELISA

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.