

**HDBS0213**

---

## Product Information

<b>Product SKU:</b>	HDBS0213	<b>Host Species:</b>	Humanized	<b>Size:</b>	100 µg
<b>Target:</b>	IL31RA	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	Q8NI17
<b>Isotype:</b>	IgG2(C14S,R16K,E20G,S21G,C102S,H147Q,R234Q,S257A,Q298E)		

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CRL; GPL; CRL3; GLMR; GLM-R; PLCA2; hGLM-R; IL-31RA; PRO21384; zcytoR17

**Common Name:** CIM 331,CIM331

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.