

HDBS0228

Product Information

Product SKU:	HDBS0228	Host Species:	Chimeric	Size:	100 µg
Target:	CD138	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P18827
Isotype:	IgG4		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: SDC1; Syndecan-1; CD138; SYND1; SDC

Common Name: BT-062, BT062, nBT062-DM4

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.