

### Product Information

<b>Product SKU:</b>	HDBS0245	<b>Host Species:</b>	Humanized	<b>Size:</b>	100 µg
<b>Target:</b>	CD22	<b>Reactivity:</b>	Human		

### Additional Information

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	P20273
<b>Isotype:</b>	IgG1		

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD22;SIGLEC2;BL-CAM;SIGLEC-2;Siglec2;SIGLEC2FLJ22814

**Common Name:** ACD22-VCMMAE, DCDT2980S, FCU2803, RG-7593, RO5541072-000

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.