

# Biotinylated Anti-CEACAM5 (Labetuzumab biosimilar) mAb

HDBS0252



## Product Information

**Product SKU:** HDBS0252      **Host Species:** Humanized      **Size:** 100 µg  
**Target:** CEACAM5      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P06731  
**Isotype:** IgG1

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CEACAM-5;CD66e;CEA;Meconium antigen 100

**Common Name:** IMMU-130, hMN-14-SN-38, hMN-14-SN-38 ADC, hMN14-CL-SN-38, hMN14-SN-38

**Recommended Dilutions:** ELISA 1:5000-10000

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.