

# Biotinylated Anti-CD40 Ligand (Ruplizumab biosimilar) mAb

HDBS0256



## Product Information

**Product SKU:** HDBS0256      **Host Species:** Humanized      **Size:** 100 µg  
**Target:** CD40 Ligand      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P29965  
**Isotype:** IgG1

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD40LG;CD154;CD40L;HIGM1;IGM;IMD3;T-BAM;TNFSF5;TRAP;gp39

**Common Name:** BG-9588, hu5c8

**Recommended Dilutions:** ELISA 1:5000-10000; Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.