

Product Information

Product SKU:	HDBS0273	Host Species:	Humanized	Size:	100 µg
Target:	ITGA4	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P13612
Isotype:	IgG4		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	CD49D;IA4
Common Name:	AN100226,AN100226,BG00002
Recommended Dilutions:	Flow Cyt 1:10
Applications:	Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.