

Biotinylated Anti-CLDN18.2 (Zolbetuximab biosimilar) mAb

HDBS0274



Product Information

Product SKU:	HDBS0274	Host Species:	Chimeric	Size:	100 µg
Target:	CLDN18.2	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	P56856
Isotype:	IgG1		

Protein Information

Background:	Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.
Synonyms:	GC-182;MAB-362;MAB362;laudiximab
Common Name:	GC-182,IMAB-362,IMAB362,claudiximab
Recommended Dilutions:	Flow Cyt 1:100
Applications:	Flow Cyt
Usage:	Research use only
Formulation:	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Storage:	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.