

HDBS0281

Product Information

| | | | | | |
|---------------------|----------|----------------------|-----------|--------------|--------|
| Product SKU: | HDBS0281 | Host Species: | Humanized | Size: | 100 µg |
| Target: | DLL3 | Reactivity: | Human | | |

Additional Information

| | | | |
|-------------------|--------------|-----------------|--------|
| Conjugate: | Biotinylated | Uniprot: | Q9NYJ7 |
| Isotype: | IgG1 | | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: SCDO1

Common Name: SC0001-SCX,SC0002,sc0001-SCX

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.