

HDBS0302

---

## Product Information

**Product SKU:** HDBS0302      **Host Species:** Homo sapiens      **Size:** 100 µg  
**Target:** CXCR4      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P61073  
**Isotype:** IgG4

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** CD184;D2S201E;FB22;HM89;HSY3RR;LAP-3;LAP3;LCR1;LESTR;NPY3R;NPYR;NPYRL;NPYY3R;WHIM;WHIMS

**Common Name:** BMS-936564, MDX-1338

**Recommended Dilutions:** Flow Cyt 1:100

**Applications:** Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.