

### Product Information

|                     |          |                      |              |              |        |
|---------------------|----------|----------------------|--------------|--------------|--------|
| <b>Product SKU:</b> | HDBS0312 | <b>Host Species:</b> | Homo sapiens | <b>Size:</b> | 100 µg |
| <b>Target:</b>      | C5AR1    | <b>Reactivity:</b>   | Human        |              |        |

### Additional Information

|                   |              |                 |        |
|-------------------|--------------|-----------------|--------|
| <b>Conjugate:</b> | Biotinylated | <b>Uniprot:</b> | P21730 |
| <b>Isotype:</b>   | IgG1         |                 |        |

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** C5a-R;C5aR;CD88

**Common Name:** IPH-5401, NN-8210, NNC-0215-0384

**Recommended Dilutions:** Flow Cyt 1:100

**Applications:** Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.