

Product Information

Product SKU:	HDBS0345	Host Species:	Humanized	Size:	100 µg
Target:	PDL1	Reactivity:	Human		

Additional Information

Conjugate:	Biotinylated	Uniprot:	Q9NZQ7
Isotype:	IgG1(K97R, L117F, L118E, P214S)		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD-L1;CD274;B7-H1;PDCD1L1;PDCD1LG1

Common Name: MEDI4736

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.