

### Product Information

**Product SKU:** HDBS0361      **Host Species:** Humanized      **Size:** 100 µg  
**Target:** ALPP      **Reactivity:** Human

### Additional Information

**Conjugate:** Biotinylated      **Uniprot:** P05187  
**Isotype:** IgG1

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** PLAP-1;PLAP

**Common Name:** A06

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.