

Biotinylated Anti-IGF-1R (Teprotumumab biosimilar) mAb

HDBS0369



Product Information

Product SKU: HDBS0369 **Host Species:** Homo sapiens **Size:** 100 µg
Target: IGF-1R **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** P08069
Isotype: IgG1(E239D, M241L)

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD221, IGFIR, IGFR, JTK13

Common Name: RO4858696-000,teprotumumab-trbw

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).Lyophilized antibodies are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.