

# Biotinylated Anti-GPNMB (Glembatumumab biosimilar) mAb

HDBS0387



## Product Information

**Product SKU:** HDBS0387      **Host Species:** Homo sapiens      **Size:** 100 µg  
**Target:** GPNMB      **Reactivity:** Human

## Additional Information

**Conjugate:** Biotinylated      **Uniprot:** Q14956  
**Isotype:** IgG2(S257A)

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** HGFIN, NMB

**Common Name:** CDX-011 (DOX), CR011, glembatumumab vedotin (ADC)

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: [techsupport@assaygenie.com](mailto:techsupport@assaygenie.com) | Web: [www.assaygenie.com](http://www.assaygenie.com)

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.