

HDBS0390

Product Information

Product SKU: HDBS0390 **Host Species:** Chimeric Humanized **Size:** 100 µg
Target: LRRC15 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** Q8TF66
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: LIB, LRC15

Common Name: PR-1498487

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.