

### Product Information

<b>Product SKU:</b>	HDBS0395	<b>Host Species:</b>	Chimeric Rabbit	<b>Size:</b>	100 µg
<b>Target:</b>	PD-1	<b>Reactivity:</b>	Canine		

### Additional Information

<b>Conjugate:</b>	Biotinylated	<b>Uniprot:</b>	NP_001301026.1
<b>Isotype:</b>	IgG		

### Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** PD1, PDCD1, CD279, SLEB2, hPD-1, hPD-I, hSLE1

**Common Name:** 4F12

**Recommended Dilutions:** ELISA 1:5000-10000, Flow Cyt 1:100

**Applications:** ELISA, Flow Cyt

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.