

Biotinylated Anti-TNFSF15 (PF-06480605 biosimilar) mAb

HDBS0417



Product Information

Product SKU: HDBS0417 **Host Species:** Humanized **Size:** 100 µg
Target: TNFSF15 **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** O95150
Isotype: IgG1(L117A,L118A,G120A)

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: TL1;VEGI

Common Name: PF 6480605、PF-6480605、RVT-3101,1D1-1.31

Recommended Dilutions: ELISA 1:5000-10000, Flow Cyt 1:100

Applications: ELISA, Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | Web: www.assaygenie.com

Copyright © 2024 Assay Genie Ltd, All Rights Reserved. All information / detail is correct at time of going to print.