

HDBS0420

Product Information

Product SKU: HDBS0420 **Host Species:** Humanized **Size:** 100 µg
Target: CRTAM **Reactivity:** Human

Additional Information

Conjugate: Biotinylated **Uniprot:** O95727
Isotype: IgG1 (K97R,L117A,L118A,E239D,M241L)

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD355

Common Name: 5A11

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.