

HDBS0436

Product Information

Product SKU: HDBS0436	Host Species: Homo sapiens	Size: 100 µg
Target: ILDR2	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: Q71H61
Isotype: IgG2(S257A)	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: C1orf32; Angulin-3

Common Name: BAY-1905254

Recommended Dilutions: Flow Cyt 1:100

Applications: Flow Cyt

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.