

HDBS0444

Product Information

Product SKU:	HDBS0444	Host Species:	Homo sapiens	Size:	100 µg
Target:	IL15	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P40933
Isotype:	IgG1(K97R, E239D, M241L)		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL-15

Common Name: HuMax-IL15, AMG 714, AMG-714, HuMax-IL-15, MAb146B7, PRV-015

Recommended Dilutions: ELISA 1:5000-10000

Applications: ELISA

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized antibodies are shipped at ambient temperature.