

HDBS0477

Product Information

Product SKU:	HDBS0477	Host Species:	Whole mAb	Size:	100 µg
Target:	TNFRSF10B	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	O14763
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: DR5;CD262;KILLER;TRICK2;TRICKB;ZTNFR9;TRAILR2;TRICK2A;TRICK2B;TRAIL-R2;KILLER/DR5

Common Name: GEN1029;GEN-1029;Hx-DR5-01;Hx-DR5-05;iDD 004;iDD-004;iDD004

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.