

HDBS0510

---

## Product Information

**Product SKU:** HDBS0510      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** C5      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** P01031  
**Isotype:** IgG1

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** C5D;C5a;C5b;ECLZB;CPAMD4

**Common Name:** CH 7092230;CH7092230;CH-7092230;RG 6107;RG-6107;RG6107;RO 7092230;RO-7112689;RO7112689;SKY 59;SKY-59;SKY59

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.