

HDBS0536

Product Information

Product SKU:	HDBS0536	Host Species:	Whole mAb	Size:	100 µg
Target:	CSF1R	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P07333
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: FMS;CSFR;FIM2;GPSC;HDLS;C-FMS;CD115;HDLS1;CSF-1R;BANDDOS;M-CSF-R

Common Name: RG-7155;RG7155;RO-5509554;RO5509554

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.