

HDBS0548

---

## Product Information

<b>Product SKU:</b> HDBS0548	<b>Host Species:</b> Whole mAb	<b>Size:</b> 100 µg
<b>Target:</b> PCSK9	<b>Reactivity:</b> Human	

---

## Additional Information

<b>Conjugate:</b> Unconjugated	<b>Uniprot:</b> Q8NBP7
<b>Isotype:</b> IgG2	

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** FH3;PC9;FHCL3;NARC1;LDLCQ1;NARC-1;HCHOLA3

**Common Name:** AMG-145;AMG145;Repatha

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.