

HDBS0550

Product Information

Product SKU: HDBS0550 **Host Species:** Whole mAb **Size:** 100 µg
Target: SNCA **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P37840
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: PD1;NACP;PARK1;PARK4

Common Name: PR-1808164;PR1808164;ABBV-0805;ABBV0805;ABBV805;ABBV-805;BAN-0805;BAN0805;BAN805;BAN-805;PD-0805;PD0805;PD-805;PD805

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.