

HDBS0557

Product Information

Product SKU: HDBS0557 **Host Species:** Whole mAb **Size:** 100 µg
Target: CD38 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P28907
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD38;ADPRC1;cADPR1;ADPRC 1;T10;cADPr hydrolase 1

Common Name: MOR 202;MOR202;MOR-202;MOR-03087;MOR 03087;MOR03087;MOR 3087;MOR-3087;MOR3087;TJ-202;TJ202;TJ 202

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.