

HDBS0573

Product Information

Product SKU: HDBS0573	Host Species: Whole mAb	Size: 100 µg
Target: CTLA-4	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: P16410
Isotype: IgG1	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD;GSE;GRD4;ALPS5;CD152;CTLA-4;IDDM12;CELIAC3

Common Name: BNT-316;BNT316;BNT 316;ONC-392;ONC392;ONC 392

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.