

HDBS0582

Product Information

Product SKU:	HDBS0582	Host Species:	Whole mAb	Size:	100 µg
Target:	IL2RA	Reactivity:	Human		

Additional Information

Conjugate:	Unconjugated	Uniprot:	P01589
Isotype:	IgG1		

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: p55;CD25;IL2R;IMD41;TCGFR;IDDM10;IL2-RA;IL-2-RA

Common Name: AU-007;AU007;AU 007;BD-8;BD8;BD 8;BDG8;BD-G8;BD G8

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.