

HDBS0590

Product Information

Product SKU: HDBS0590	Host Species: Whole mAb	Size: 100 µg
Target: IL17A	Reactivity: Human	

Additional Information

Conjugate: Unconjugated	Uniprot: Q16552
Isotype: IgG4	

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: IL17;CTLA8;IL-17;ILA17;CTLA-8;IL-17A

Common Name: LY-2439821;LY2439821;LY 2439821;Taltz

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.