

HDBS0591

Product Information

Product SKU: HDBS0591 **Host Species:** Whole mAb **Size:** 100 µg
Target: C5AR1 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P21730
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: C5A;C5AR;C5R1;CD88

Common Name: HIB210;HIB-210;HIB 210;MOR044254;MOR-044254;MOR 044254;MOR210;MOR-210;MOR 210;TJ 210;TJ-210;TJ210;TJ210001;WBP2191;WBP-2191;WBP 2191

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.