

HDBS0595

Product Information

Product SKU: HDBS0595 **Host Species:** Whole mAb **Size:** 100 µg
Target: IL13 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P35225
Isotype: IgG4

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: P600;IL-13;NC30

Common Name: DRM 06;DRM-06;DRM06;DRM 6;DRM-6;DRM6;LY 3650150;LY-3650150;LY3650150;MILR1444A;MILR-1444A;MILR 1444;ARG 3637;ARG-3637;ARG3637;RO 5490255;RO-5490255;RO5490255;TNX 650;TNX-650;TXN650

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.