

HDBS0602

Product Information

| | | |
|------------------------------|--------------------------------|---------------------|
| Product SKU: HDBS0602 | Host Species: Whole mAb | Size: 100 µg |
| Target: TGFB1 | Reactivity: Human | |

Additional Information

| | |
|--------------------------------|------------------------|
| Conjugate: Unconjugated | Uniprot: P01137 |
| Isotype: IgG4 | |

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CED;LAP;DPD1;TGFB;IBDIMDE;TGFbeta;TGF-beta1

Common Name: SRK-181;SRK 181;SKR181

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.