

HDBS0611

---

## Product Information

**Product SKU:** HDBS0611      **Host Species:** Whole mAb      **Size:** 100 µg  
**Target:** Her2      **Reactivity:** Human

---

## Additional Information

**Conjugate:** Unconjugated      **Uniprot:** P04626  
**Isotype:** IgG1

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** NEU;NGL;HER2;TKR1;CD340;HER-2;VSCN2;MLN 19;MLN-19;c-ERB2;c-ERB-2;HER-2/neu;p185(erbB2)

**Common Name:** PF-06741086;PF06741086;PF-6741086;PF6741086

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.