

HDBS0620

Product Information

Product SKU: HDBS0620 **Host Species:** Whole mAb **Size:** 100 µg
Target: CD40 **Reactivity:** Human

Additional Information

Conjugate: Unconjugated **Uniprot:** P25942
Isotype: IgG1

Protein Information

Background: Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

Synonyms: CD40;p50;Bp50;CDW40;TNFRSF5

Common Name: ADC 1013;ADC-1013;ADC1013;JNJ 7107;JNJ7107;JNJ-7107;JNJ-64457107;JNJ64457107;JNJ 64457107;Vanalimab

Recommended Dilutions: -

Applications: -

Usage: Research use only

Formulation: Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

Storage: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.