

HDBS0623

---

## Product Information

<b>Product SKU:</b>	HDBS0623	<b>Host Species:</b>	Whole mAb	<b>Size:</b>	100 µg
<b>Target:</b>	CD73	<b>Reactivity:</b>	Human		

---

## Additional Information

<b>Conjugate:</b>	Unconjugated	<b>Uniprot:</b>	P21589
<b>Isotype:</b>	IgG1		

---

## Protein Information

**Background:** Research grade biosimilar. Not for use in therapeutic or diagnostic procedures for humans or animals.

**Synonyms:** NT;eN;NT5;NTE;eNT;CD73;E5NT;CALJA

**Common Name:** COR-004;COR004;COR 004;COR-4;COR4;COR 4;CPI-006;CPI006;CPI 006;CPI-6;CPI6;CPI 6;CPX-006;CPX006;CPX 006;CPX-6;CPX6;CPX 6

**Recommended Dilutions:** -

**Applications:** -

**Usage:** Research use only

**Formulation:** Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.

**Storage:** Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.